

Estd. 1966

®

DESCO

DELIVERING TO THE WORLD

Autoloader Collapsible Stretcher (Ambulance Stretcher) (MODEL NO : STAS 101)

Make in India



Decide with Confidence

Estd. 1966

®

DESCO

DELIVERING TO THE WORLD

Hospital Autoloader Collapsible Stretcher (Ambulance Stretcher)



www.medicalhospitalfurnitures.com

www.descoonline.in

www.descomedicalindia.com

Autoloader Collapsible Stretcher (Ambulance Stretcher) Supplied by Medical Hospital Furniture

Specs :-

* This automatic loading stretcher is our latest product, it's made of high-strength aluminum alloy thick-wall pipes.

* The bed surface is made of foamed cushion and the adjustable back making the wounded feel comfortable.

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)

Estd. 1966



DESCO

DELIVERING TO THE WORLD

Autoloader Collapsible Stretcher (Ambulance Stretcher) (MODEL NO : STAS 101)



- * The mechanical foldaway legs are controlled by the handles on both ends of the stretcher.
- * And the legs can be folded under the stretcher automatically when loading into the vehicle.
- * Also the stretcher can be locked by the fixing devices when on the ambulance car.
- * This type of stretcher is characterized by its being light.
- * Durable antiseptic easy for sterilization.

Technical Parameter :

- * High Position (L x W x H) : 1900 x 550 x 920mm.
- * Low Position : 1900 x 550 x 250mm.
- * Packing Size (1 pc) : 1950 x 620 x 255mm
- * Max Angle of the Back : 75°
- * Load Bearing : ≤ 159 kg
- * N.W. (1 pc) : 40 kg
- * G.W. : 46 kg

We are Exporters and Suppliers of Autoloader Collapsible Stretcher (Ambulance Stretcher) at competitive prices from India.

MODEL NO:

STAS 101 Aluminium Anodized