

Estd. 1966

®

DESCO

DELIVERING TO THE WORLD

Electric Luxury Obstetric Table (MODEL NO : DTEL 105)

Make in India

CE



Decide with Confidence

Estd. 1966

®



DELIVERING TO THE WORLD

Hospital Electric Luxury Obstetric Table



www.medicalhospitalfurnitures.com

www.descoonline.in

www.descomedicalindia.com

Electric Luxury Obstetric Table Manufactured by Medical Hospital Furniture

Specs :-

This Table is Suitable for Examination Diagnosis & Operation for Gynaecology, Urological Surgery.

Operating Positions

Three Sectional upholstered tabletop mounted on strong S.S. tubular frame ensures maximum comfort to the patient.

- Maintenance free lifting system with heavy duty ball bearings provide silent & jerk free movements.
- Built-in battery backup for emergency use only.
- Fully Electric Table with Operating Positions.

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)

Estd. 1966



DESCO

DELIVERING TO THE WORLD

Electric Luxury Obstetric Table (MODEL NO : DTEL 105)



Height Adjustment

- Trendelenburg & Reverse Trendelenburg.
- Back Section Lift.
- SS 304 Grade steel Used.
- Leg section is easily detachable and floor locking is manual.
- Table is suitable for examination diagnosis and operation for Gynecology, Urological Surgery.
- Actuators used Ti Motion (Taiwan).
- Dual Access control-Hand held Remote & Foot Switch.
- Specially design base with 'U' Cut providing extra leg room while examination.
- * Table Size (LXW) : 191 x 62 cm.
- * Min. Height – 78 cm.
- * Max. Height – 102cm.
- * Trendelenburg / Rev. Trendelenburg – 25 Degree.
- * Back Section - 70 Degree.

Standard Accessories

- * 'L' Shaped Anaesthetic Screen : 1 Pc.
- * Padded Shoulder Support : 1 Pair.
- * Radio-translucent Armboard with pad : 2 Pcs.
- * Goepel Knee Crutch with straps : 1 Pair.
- * PU Foam Mattress : 1 Set.
- * Water proof table top cover : 1 Pc.
- * S.S. Douche Tray with drain pipe : 1 Pc.

We are Manufacturer, Exporters and Suppliers of Electric Luxury Obstetric Table at competitive prices from India.

MODEL NO:

DTEL 105 Electric Luxury Obstetric Table