

Estd. 1966

®

DESCO

DELIVERING TO THE WORLD

Electric Obstetric Table (MODEL NO : DTEL 101)

Make in India

CE



Decide with Confidence


Estd. 1966

®

DESCO

DELIVERING TO THE WORLD

Hospital Electric Obstetric Table



www.medicalhospitalfurnitures.com

www.descoonline.in

www.descomedicalindia.com

Electric Obstetric Table Manufactured by Medical Hospital Furniture

Specs :-

Obstetrics operating room table perfect for all kinds of gynecological operations.

The cover for base and column is made of 304 stainless steel, good performance on anti-rust and easy to clean.

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)

Estd. 1966



DESCO

DELIVERING TO THE WORLD

Electric Obstetric Table (MODEL NO : DTEL 101)



Technical Data

- Length of Table : 1600 mm.
- Width of Table Top : 550 mm.
- Minimum Height (without Mattress) : Min. 750 mm.
- Maximum Height (without Mattress) : Max. 1000 mm.
- Lateral Tilt : $\pm 20^\circ$.
- Back Section : $\pm 70^\circ$.
- Floor Locking System : Manual.

Features

- Functional and versatile.
- Flexibility and precision.
- Hygienic design.
- Easy cleaning.
- User- friendly.
- STEEL QUALITY 304 MEDICAL GRADE STAINLESS STEEL.

Standard Accessories

- * 'L' Shaped Anaesthetic Screen : 1 Pc.
- * Padded Shoulder Support : 1 Pair.
- * Radio-translucent Armboard with pad : 2 Pcs.
- * Goepel Knee Crutch with straps : 1 Pair.
- * PU Foam Mattress : 1 Set.
- * Water proof table top cover : 1 Pc.
- * S.S. Douche Tray with drain pipe : 1 Pc.

We are Manufacturer, Exporters and Suppliers of Electric Obstetric Table at competitive prices from India.

MODEL NO:

DTEL 101 Electric Obstetric Table