

Estd. 1966

DESCO

DELIVERING TO THE WORLD

Obstetric Table Deluxe (MODEL NO : DTEL 110)

Make in India



Decide with Confidence

Hospital Obstetric Table Deluxe

Estd. 1966



DELIVERING TO THE WORLD



www.medicalhospitalfurnitures.com

www.descoonline.in

www.descomedicalindia.com

Electric Luxury Obstetric Table Manufactured by Medical Hospital Furniture

Specs :-

This Table is specially designed for Multi Function Electric Obstetric Table with Three electric movements suitable for Gynaecological diagnosis & examination , Child Birth & Ceasarian Operation

Operating Positions :-

* Height adjustment, Trendelenburg, Reverse Trendelenburg & back section is smoothly controlled by handheld remote as well as foot switch using linear motors.

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)

Estd. 1966



DESCO

DELIVERING TO THE WORLD

Obstetric Table Deluxe (MODEL NO : DTEL 110)



- * 3 section upholstered table top ensures maximum comfort to the patients.
- * Leg section can be completely slide in under the middle section for lithotomy position or slide out into full length table as per delivery requirements.
- * Maintenance free lifting system with hi- precision bearings ensuring highest rigidity & smooth up & down movements.
- * Complete with SS 304 Grade Stainless Steel and all standard accessories.

Technical Specifications

- Table Top Size (LXW) : 1320 x 660 mm
- Leg Section (LXW) : 550 x 600 mm
- Table Top Size with Leg Section : 1870 x 660 mm
- Min Height : 78 cm, Max Height : 1030 mm
- Trendelenburg / Reverse Trendelenburg : 15 Degree/15 Degree
- Back Section Adjustment : +0 Degree/60 Degree
- Load capacity : 180 Kg.
- Power supply : 220 VAC -50Hz
- Drives : Linear Actuators
- Battery Backup : Upto 7 Cycles (Emergency use only)

Standard Accessories

- * 'L' Shaped Anaesthetic Screen : 1 Pc.
- * Padded Shoulder Support : 1 Pair
- * Radio-translucent Armboard with pad : 2
- * Goepel Knee Crutch with straps : 1 Pair
- * PU Foam Mattress : 1 Set
- * Water proof table top cover : 1 Pc.
- * S.S. Douche Tray with drain pipe : 1 Pc.

We are Manufacturer, Exporters and Suppliers of Obstetric Table Deluxe at competitive prices from India.

MODEL NO:

- DTEL 110** Obstetric Table Deluxe Hydraulic
- DTEL 111** Obstetric Table Deluxe Semi Electric
- DTEL 112** Obstetric Table Deluxe Electric